


Validation Master Plan for Sterilising-Grade Filters



Pharmaceutical Application and
Validation Services



pure solutions



Factors Influencing Sterility Assurance Level

- Filter performance
- Filtration process design
 - ◆ System sizing and engineering
 - ◆ Standard Operating Procedure
- Validation of the filtration process
- Training of operators
- Process control and monitoring





Perspectives

MILLIPORE

focuses on increasing Sterility Assurance Level

- ◆ Reliable sterilizing-grade filters
- ◆ Rigorous QA/QC system
- ◆ Application expertise and validation services

Drug manufacturer

focuses on Potency Efficacy Purity Safety Compliance

- ◆ High quality filtration products
- ◆ Robust filtration process
- ◆ Sterility assurance

MILLIPORE



Sterilising-grade filter

Definition

➤ FDA Guideline on Sterile Drug Products Produced by Aseptic Processing 1987

- ◆ A sterilizing filter is one which, when challenged with the microorganism *Brevundimonas diminuta*, at a minimum concentration of 10^7 organisms per cm^2 of filter surface, will produce a sterile effluent.
- ◆ Correlate filter performance with filter integrity testing.

2004

- ◆ A filter which, when appropriately validated, will remove all microorganisms from a fluid stream, producing sterile effluent.
- ◆ A product filter's integrity test specification should be consistent with data generated during bacterial retention validation studies.

➤ PDA technical report 26 1998

- ◆ A physical integrity test is meaningful only when it can be related to specific filter retention characteristics



Sterilising-grade filter

How does a membrane retain ?

Retention Mechanisms

Process Fluid





Size Exclusion

Adsorption

Sieving

Entrapment



MILLIPORE




Sterilising-grade filter


Single layer & controlled pore size

Size exclusion

controlled pore size distribution
less dependant on process conditions





Surface view



Cross section
view

Adsorption

added thickness
double layer structure



MILLIPORE



Validation Master Plan

- Validation of filter performance
 - ◆ Qualification of sterilising-grade filter
 - ◆ Specific validation study
- Validation of filtration process
 - ◆ Filter sterilisation
 - ◆ Filter integrity testing
 - ◆ Process control and monitoring
- Training and certification
- Quality Control & Quality Assurance system



MILLIPORE


MILLIPORE Manufacturing & Quality Assurance System



pure solutions



Membrane Casting





MILLIPORE



Membrane Quality Assurance


- Bacterial retention
- Bubble point
- Thickness
- Flow time
- Pyrogenicity
- Toxicity



MILLIPORE



Lamination & Pleating



POLYPROPYLENE
BAND

HEATING BLOCKS

MEMBRANE

POLYPROPYLENE


MEMBRANE

POLYPROPYLENE

PLEATING

PACK


CUTTER



MILLIPORE





Seaming & End capping



MILLIPORE



Serial numbering for complete traceability & Stacking



MILLIPORE



Device Quality Assurance

- Bacterial retention
- Bubble point
- 100% device integrity testing
- Liquid flow rate
- Extractables
- USP bacterial endotoxin
- Mechanical stress (thermal, hydraulic)
- Steaming resistance



MILLIPORE



Position Statement

Millipore ensures sterilising-grade performance for filters by using a three step integrity testing strategy:

1. All membranes used in Millipore sterilising-grade filters are BP tested
2. After filter devices are fabricated, Millipore performs a proprietary integrity test on all products
3. Finally, Millipore provides customers with integrity test specification

MILLIPORE



100 % In-process Device Integrity Testing

Integral membrane (*wetted*)

Damaged membrane (*wetted*)

Air
diffusion

spec = 10 ml/min

ratio = 1.8

5 ml/min

4 ml/min Defect

9 ml/min

Low Solubility Gas
diffusion

spec = 0.2 ml/min

ratio = 41

0.1 ml/min

4 ml/min Defect

4.1 ml/min



Requirements for increased SAL

➤ PDA TR n° 26

- ◆ “Integrity testing alone is insufficient to assure the sterility of the filtrate. At least two other elements must be in place:
- ◆ the production controls and quality assurance systems used by the filter manufacturer...,
- ◆ and the validation studies used to show that a particular combination of product, processing conditions and sterilizing grade filter will meet the requirements of the bacterial challenge test”.

MILLIPORE


Specific Validation Study



pure solutions




pure solutions



MILLIPORE




Validation of Filter Performance





There is NO absolute filter

Factors that can influence filter retention





Bacteria retention testing

Regulatory requirements

- FDA Guideline on Sterile Drug Products Produced by Aseptic Processing 1987, 2002
 - ◆ Validation should include microbiological challenges
 - ◆ Acceptable challenge is at least 10^7 cfu/cm²
 - ◆ The micro-organism should simulate the smallest micro-organism that may occur in production
 - ◆ Laboratory experiments should simulate production conditions: flow rates, pressures, temperature, time



Bacteria retention testing

cGMP advice

PDA technical report 26

- Filter membrane lots used for bacterial retention should have a pre-filtration water wet physical integrity test value at or near the filter manufacturers production limit
- The integrity test will help establish the similarity of the test filters to the filters validated to retain a bacterial challenge under process-related conditions



Extractables testing

Regulatory requirements

- FDA Guideline for Submitting Documentation for Sterilization Process Validation
 - ◆ Any effects of the filter on the product formulation should be described (adsorption or extractables)
- FDA Human Drug cGMP Notes 9/94
 - ◆ Extractables cannot be detected because drug product interferes with methods and quantities are very low
 - ◆ Use appropriate methods and solvents to obtain the amount of extractables per filter.
 - ◆ Show the identity, quantity, and toxicity of the extractables.



Filter integrity testing


- Water is the standard wetting fluid
 - ◆ Drug component binding reduce surface tension
 - ◆ Post-use test failures

- Testing with the drug product as the wetting agent
 - ◆ Facilitates pre-use post-SIP testing
 - ◆ Avoids rinsing the membrane after filtration
 - ◆ Reduces the risk of post-use integrity testing failures




Compatibility & Adsorption



- Expose filters to the drug solution and the processing conditions
- Evaluate adsorption
 - ◆ Analytical assay for final product release
- Check compatibility
 - ◆ Confirm filter integrity
 - ◆ Also confirmed during
 - Microbial retention
 - Extractables testing



Validation of the Aseptic Process




pure solutions





pure solutions




MILLIPORE



Design of filtration process

- Design of filtration system
 - ◆ P&ID and flow schematics
 - ◆ Functional description
- Development of SOPs
 - ◆ Filter SIP & autoclaving
 - ◆ Filter integrity testing
- Installation qualification
- Operational qualification
- Performance qualification



MILLIPORE



Validation of filter integrity tester

Regulatory requirements

- EU GMP
 - ◆ Equipment which are critical for the quality of the products shall be subjected to appropriate qualification
- FDA Guideline on general principles of process validation
 - ◆ The manufacturer shall validate computerized systems that monitor and/or control the manufacturing process



MILLIPORE



Validation of filter sterilisation

Regulatory requirements

➤ FDA Aseptic guideline 1987

- ◆ In-line filters cause a significant pressure differential and temperature drop
- ◆ Use suitable biological indicators at appropriate downstream locations
- ◆ Filter SIP validation should also include measurements of temperature and pressure

➤ EU GMP Annex 1 1997

- ◆ ... Any sterilisation process ... should be demonstrated by thermometric means and by biological indicators where appropriate



MILLIPORE



Bioburden monitoring

Regulatory requirements

- FDA Aseptic guideline 1987
 - ◆ Assure that actual bioburden does not contain micro-organisms of a size and/or concentration that would reduce sterility assurance
- FDA Recommendations for submitting documentation for sterilization process validation
 - ◆ Describe the program for routinely monitoring bioburden
 - ◆ Give information concerning the number, type, and resistance of bioburden organisms



Bioburden monitoring

EU Regulatory requirements

- Manufacture of the finished dosage form
EU Guideline 3AQ2a
 - ◆ For sterilisation by filtration the maximum acceptable bioburden...in most situations, 10 cfu's/100 ml
 - ◆ If this requirement is not met, it is necessary to use a pre-filtration through a bacteria-retaining filter to obtain a sufficiently low bioburden



Validation of bioburden monitoring

➤ Equipment / device / media

Qualification

Environmental

Product bioburden

➤ Method validation

Personnel

➤ Monitoring program

(Sterility testing)

◆ Absence of bioburden

➤ Maximum microbial limits



Process controls

- FDA Aseptic guideline 1987, 2002
 - ◆ The total time for product filtration ... should be limited to an established maximum

- EU GMP Manufacture of sterile medicinal products
Annex 1 1997
 - ◆ The time taken to filter a known volume of bulk solution and the pressure difference to be used across the filter should be determined during validation and any significant differences from this should be noted and investigated



Training

➤ EU GMP

- ◆ Provide training for personnel whose activities could affect the quality of the product
- ◆ Personnel should receive training appropriate to the duties assigned to them

➤ US GMP

- ◆ Each responsible person shall have the training and experience to perform the assigned functions





Validation documentation

VMP




Validation Rationale



Access Questionnaire



Protocol



Report



Validation Master file

MILLIPORE



Sterile filtration validation Responsibilities

Test type	Millipore access services	Drug manufacturer
Microbial retention testing Extractables Compatibility Product-specific integrity-test data	Protocol Laboratory execution Report	Approve protocols Maintain documentation Process controls
Adsorption	Protocol	Execution, report
Integrity testers IQ/OQ	Protocol Field execution & report	Approve protocols Maintain documentation
Validation of integrity-testing PQ, SOP, training	Operator training and certification	PQ protocol Field execution, report
Validation of filter-sterilisation PQ, SOP, training	Protocol Operator training	Field execution, report Documentation




Bioburden monitoring validation

Responsibilities

Test type	Millipore access services	Drug manufacturer
Method development Product bioburden Sterility testing	Protocol Laboratory execution Report	Approve protocol Maintain documentation
Method validation Equipment IQ/OQ Microbiological procedures Training	Protocol Lab execution & report (Documentation) Certification	Approve protocols (Lab execution & report) Maintain documentation




Qualification route





pure solutions



MILLIPORE



References

EMEA

eudraportal.eudra.org
www.emea.eu.int

FDA

Aseptic guideline
HcGMP notes

www.fda.gov

US GMP

21CFR212

www.access.gpo.gov

EU GMP

Eudralex volume 4

<http://pharmacos.eudra.org/F2/home.html>

PDA

Technical report 26

www.pda.org